H.I.T.[®] - Hook Inserting Tack System


H.I.T.[®]
Insertion
Pliers

H.I.T.[®]
Hooks

H.I.T.[®]
Punch
Pliers

 **Manufacturer:**
Dentsply International Raintree Essix
7290 26th Court East
Sarasota, FL 34243 USA
Tel: 800.883.8733
www.essix.com

Emergo Europe
Prinsessegracht 20
2514 AP, The Hague
The Netherlands

R_x
ONLY

CE
0120

H.I.T.[®] – Hook Inserting Tack System

ENG

Function:

The H.I.T.[®] – Hook Inserting Tack System is designed to insert external hooks onto thermoplastic aligners for the purposes of attaching elastics.

Contraindications:

The H.I.T.[®] – Hook Inserting Tack System is not designed to place force points for purposes of tooth movement.

⚠ Warnings:

Treatment range using the H.I.T.[®] – Hook Inserting Tack System is not to exceed 1mm of movement per month per tooth. Excessive rate of movement may cause root resorption.

Precautions:


None known.

Adverse Reactions:

None known.

Step by Step Instructions:

1. Pliers must be sterilized before each use. See sterilization guidelines below.
2. Space must be created in the aligner before placing the hook. Use the Thermoplier[®] Large Circle Pliers (Gray) to raise the plastic away from the enamel of target tooth.
3. Using the H.I.T.[®] Punch Pliers, punch a hole in the center of the raised surface of the aligner.
4. Place the thinner, tapered side of the H.I.T.[®] Hook into the head of the H.I.T.[®] Insertion Pliers. Then, insert the hook from the outside of the aligner into the hole.

H.I.T.[®] Pliers Purpose/Indication:

Item	Purpose/Indication	Position on Aligner
H.I.T. [®] Punch Pliers	To punch a hole into which the hook is inserted	In the center of raised surface on target tooth
H.I.T. [®] Insertion Pliers	To insert a hook	In hole punched with H.I.T. [®] Punch Pliers on target tooth
H.I.T. [®] Hooks	To attach elastics	In hole punched with H.I.T. [®] Punch Pliers on target tooth

Reference the MTM[®] Clear•Aligner In-Office Quick Start Guide for detailed examples of use.

Pliers Cleaning and Sterilization Guidelines:

Clean, rinse and dry before sterilization.

Sterilize by autoclave/moist heat at 121°C / 250°F for 30 minutes (gravity displacement cycle). Minimum recommended drying time is 30 minutes.

Pliers Maintenance Suggestions:

AFTER sterilizing, lubricate SPARINGLY ONCE A WEEK with dental handpiece oil. One drop or light spray is sufficient; wipe off excess lubricant. The major cause of staining or discoloring is OVERUSE of lubricants or “sticky joint” treatments. Lubricants, waxes and “sticky joint” remedies will vaporize inside the chamber and saturate all other items placed inside.

If you are experiencing brown or yellow staining, scrub with ISOPROPYL alcohol and steel wool (type A extra fine).

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