USER MANUAL

Essix® SelectVac® Vacuum Machine

TABLE OF CONTENTS

Manual By Language

English.....2-5
 Spanish.....6-9
 French.....10-13
 Italian.....14-17
 German.....18-21

Heating Times.....22

General Information

Intended Use.....2
 Packaging Contents.....2
 Setting Up.....2
 Parts Identification.....3
 Instructions for Basic Operation.....3
 Troubleshooting.....3
 Cleaning and Maintenance.....4
 Repairs and Returns.....4
 Warranty Terms and Conditions.....4
 Technical Specifications.....5
 Exploded View.....5

INTENDED USE

This unit is intended to be used in both the dental office and laboratory to manufacture thermoplastic appliances.

PACKAGING CONTENTS

- Power Cable
- Glass Pellets
- 4mm Allen Wrench
- Cast Rack
- User Manual
- Essix® SelectVac® Machine

SETTING UP

The Essix® SelectVac® Machine requires minimal preparation for use. A few simple steps are necessary before operating the unit:

1. **Handles:** Attach the small handle to the front of the heater cover. Install the pivot brake. Then, screw the sheet support handle (the knob with the label) into the right hand side of the machine. Turn it towards you to lower the plate and start the vacuum. The blank knob screws into the left hand side and remains fixed during operation.
2. **Placement:** The Essix® SelectVac® Machine should be placed on a flat, open counter, within 4' of a standard properly grounded electrical outlet. Keep area free from moisture and excessive heat. **CAUTION:** The unit should not be placed closer than 6" from walls, cabinets or other equipment while in operation because of the temperature emitted from the heater. Under no circumstances should this unit be operated while located under or in a cabinet. Do not block bottom vent.
3. **Electric Cord:** Always ensure the machine is set up for the proper voltage (110V or 220V). Plug the female end of the cord into the electrical inlet located on the back of the unit and the male end into a properly grounded electrical outlet. The unit's power cord should not be placed near the heater or any other hot surface. When setting up the unit, do not block the inlet power cord so the unit can be disconnected if needed.


PARTS IDENTIFICATION

Figure 01

- | | |
|---------------------------------------|--------------------------------|
| 01 Heating unit handle | 07 Vacuum chamber |
| 02 Heating unit | 08 Main switch |
| 03 Sheet support handle | 09 Sheet lock handle |
| 04 Sheet support | 10 Digital timer |
| 05 Column | 11 Rotary sheet support handle |
| 06 Combined cast rack (cup/flat type) | 12 Pivot Brake |

Figure 02

- | |
|------------------------|
| 13 Rotating lock screw |
| 14 Sheet lock |
| 15 Side lock pins |

INSTRUCTIONS FOR BASIC OPERATION

1. Swing the heating unit to the rear. Raise the plastic support plate and lock in position by tightening the right-hand knob.
2. Position the model in the center of the model support plate. For taller models, flip the model support plate over and fill cup with supplied glass pellets. **NOTE:** If utilizing two-sided heating function, do not place model on base until after heating first side of plastic.
3. Place sheet of plastic on the upper support plate and secure by placing the locking ring over plastic and rotating clockwise. Swing the heater into position directly over the plastic.
4. Press the on / off button to turn on timer. Press set button once to modify timer setting, then press adjust button repeatedly to change first digit of display. Press set button to advance to next digit. Repeat until time is correct, then press set to confirm. A beep will signal that the timer is set and ready for operation.
5. Flip the heating / vacuum power switch to turn on the heater and start the timer. Observe the plastic as it heats. Follow directions for the Essix® Plastic you are using. Watch for sag depth, if plastic instructions specify. Excessive sag will result in a thin appliance. In extreme cases, the plastic may be ingested into the vacuum motor causing damage to the motor.
6. Once the timer reaches zero, it will begin beeping. To start thermoforming, twist the handle on the right side counterclockwise to allow the plastic plate to drop over the model. The vacuum will automatically start once the plastic is lowered. Allow the machine to run for the required time to cool the plastic, then flip the heating / vacuum switch off to stop the vacuum. The heater is automatically shut off once the vacuum begins.
7. Using the heating unit handle, swing the heating element to the rear of the machine. **CAUTION:** the surface of the heating unit will become hot with use and is a potential burn hazard. Contact handle only.

INSTRUCTIONS FOR TWO-SIDED HEATING

1. The machine is shipped with the rotating plate locked by a screw (Item 12, Figure 02). Prior to use of two-sided heating, loosen this screw using the supplied 4mm Allen wrench. The screw does not need to be removed.
2. Follow steps 1-4 of the basic operation process, making sure not to place the model on the plate yet.
3. Using the timer, enter the total time that will be required to heat both sides of the plastic according to the manufacturer's directions, allowing some time for flipping the plastic. The sheets may be thermoformed at any time before the timer reaches zero so it is better to have too much time than not enough.
4. Once the overall heating time has been entered, switch on the heater to begin the heating process.
5. Once the first side has been heated for the required amount of time, use the handle to swing the heating element 90° to the left.
6. Using the handle, rotate the plastic support plate **counter-clockwise** by 180° until you feel a click, thus exposing the other side of the plastic to the heater.
7. Using the handle, return the heating element into position over the plastic.
8. Heat the second side of the plastic for the required amount of time.
9. Repeat steps 5 through 7 with the exception of rotating the plastic **clockwise** to return it to the upright position.
10. If specified, continue the final heating of the plastic until it is ready for thermoforming. Once the plastic is at the ideal forming temperature, follow the steps of the basic forming process until the appliance is formed. **CAUTION:** Never lower the plastic when the plate has been rotated in the counter-clockwise direction. Always return the plastic to the original upright position before forming.

TROUBLESHOOTING

1. **Insufficient detail:** This is often caused by improper heating or vacuum time. If your adaptation does not show enough detail, allow the material to heat longer prior to forming and/or increase vacuum time. Refer to the directions included in your Essix® Brand Plastics for proper heating times.
 2. **Material creasing or tearing:** This is a result of excessive heat. Reduce the heating time. If forming several appliances in succession, unit will retain heat and require less heating time than normal. Adjust accordingly.
 3. **Material sticking to the model:** Spray the model with Trim-Rite® Release Agent prior to forming.
 4. **Material sticking to material plates:** Allow material to cool longer. Essix Freeze Spray® Coolant may be utilized to expedite this process.
 5. **Material stuck in the vacuum holes:** This can occur if the material is overheated and pulled into and through the vacuum holes. Allow materials to cool completely. Attempt to pull the material from the platform. If the material is not protruding through the vacuum holes, gently pry the material off. Excessive overheating of the material may cause permanent damage to the motor and require return for non-warranty repair.
 6. **Heater does not function:** Ensure proper power cord attachment. Ensure the power switch has been pressed to power the unit. If the heater still does not function, check the vacuum. If neither work, check the fuse located near the unit's power inlet. A spare fuse is provided with the unit. (see below for fuse replacement instructions).
- If the unit remains non-operational, please call Dentsply Sirona Orthodontics at 1-800-883-8733 for further assistance.

CLEANING AND MAINTENANCE

1. A slightly damp cloth can be used as necessary to wipe down the unit. No other cleaning should be required.
2. Check the power cord for wear or damage (ie. heat damage, fraying, etc.) at each use and examine the plug to ensure that the connections are fully seated. Never allow the heater to be placed over the power cord for extended periods as this may cause damage.
3. Keep the heating coil free from moisture and never touch with tools or metallic objects.
4. **Fuse Replacement:** The power inlet on the back of the unit incorporates a fuse. The fuse is rated for proper operation under normal conditions. If the fuse needs to be replaced, an extra fuse is included for your convenience. It is located in the fuse holder. First, remove the cord assembly from the unit by unplugging from the electrical outlet and from the unit itself. Use a flat bladed screwdriver to gently pry the fuse holder from the inlet at the area indicated by the fuse

symbol. Remove the fuse and replace it with the one located in the holder. It is important to obtain spare fuses of the same type for replacement (12 Amp fuse for the 110 volt unit and 7 Amp for the 220V unit). A different type of fuse may cause failure of the unit and will void warranty. Replacement fuses may be ordered from Dentsply Sirona Orthodontics at 1-800-883-8733.

- Only the finest quality materials are incorporated in the Essix® SelectVac® Machine. However, after years of use, the heater or some other part may need replacement. If a part requires replacement, please contact Dentsply Sirona Orthodontics at 1-800-883-8733.

REPAIRS AND RETURNS

For repair procedures, contact: Dentsply Sirona Orthodontics, 7290 26th Court East, Sarasota, FL 34243, USA.

DO NOT pack the Essix® SelectVac® Machine in styrofoam chip material for shipping. The material may lodge into the heater mechanism and become a fire hazard. Remove handles prior to packing the unit. **If you have questions or need technical assistance, please contact the Dentsply Sirona Orthodontics Customer Service Center at 800-883-8733.**

WARRANTY TERMS AND CONDITIONS

The Essix® SelectVac® Machine is warranted for its designated use, in accordance with the instructions, against original defects in material and workmanship for a period of 12 months* from date of purchase. This warranty extends only to the original purchaser. *Warranty void if unit is tampered with.

NOTE: If the equipment is not used in the manner specified by the manufacturer, the protection provided by the equipment may be impaired. Gaskets, fuses, power cord and castings are excluded from warranty guidelines.


EXPLODED VIEW

(facing page)

	Description	Qty.	110V Code	220V Code
01	Motor Housing	1	FINV1138	FINV1138
02	On/Off Switch	1	CINT0507	CINT0507
03	Convex Head Philips Screw M3.5 X 8mm	2	CPAR1004	CPAR1004
04	Washer Ø5 X 12.3 X 1.2mm	1	CARR0858	CARR0858
05	Microswitch	1	CMCR0057	CMCR0057
06	Motor Housing Lid	1	FTAP0826	FTAP0826
07	Combined Cast Rack (cup/flat type)	1	FPTM0858	FPTM0858
08	Rotary Sheet Support	1	FANE1061	FANE1061
09	Heating Unit/Sheet Support Handle	3	CCAB1345	CCAB1345
10	Sheet Lock Handle Screw	1	FPAR0919	FPAR0919
11	Sheet Lock	1	FANE1058	FANE1058
12	Pivot Handle	2	FCAB1135	FCAB1135
13	Right Pivot Handle Screw	1	FPAR1137	FPAR1137
14	Heating Unit/Rotary Sheet Support Handle Screw	2	FPAR0918	FPAR0918
15	Heating Element Housing	1	FINV1059	FINV1059
16	Convex Head Philips Screw M4 X 16mm	2	CPAR0541	CPAR0541
17	Isolator	2	CISO1256	CISO1256
18	Heating Element Reflector	1	FREF0842	FREF0842
19	Convex Head Allen Screw M5 X 16mm	2	CPAR1384	CPAR1384
20	Heating Unit Bush	1	FBCH0937	FBCH0937
21	Pivot Bush	1	FBCH0938	FBCH0938
22	Heating Element Ø10 X ØEXT 120 127V-450W	1	CRES1227	CRES1228
23	Bearing spacer	1	FESP0841	FESP0841
24	Elastic Ring Ø23.5 X 1mm	1	CANE1229	CANE1229
25	Bearing 608 ZC3 Ø22 X Ø8 X 7	1	CROL1225	CROL1225
26	Bearing Axle	1	FEIX0840	FEIX0840
27	Headless Hexagonal Screw M10X10mm	1	FPAR0868	FPAR0868
28	Spring	1	CMOL0047	CMOL0047
29	Lock Pin	1	FPIN0856	FPIN0856
30	Pivot	1	FPIV1060	FPIV1060
31	Hexagonal Column	1	FCOL0921	FCOL0921
32	End Course Pin	1	FPIN0429	FPIN0429
33	Electric Cable NEMA 5-15P	1	CCEL0498	CCEL1669
34	Universal Motor 1400W	1	CMOT1428	CMOT1429
35	Left Pivot Handle Screw	1	FPAR1136	FPAR1136
36	O'ring OR1-249	1	COR11435	COR11435
37	Lower column bush	1	CBCH1363	CBCH1363
38	Pivot Brake 15MM	1	CPIN1494	CPIN1494
39	Fuse 5 X 20mm (12A X 250V)	2	CFSV1192	CFSV1191
40	Socket with Fuse Holder	1	CTOM1261	CTOM1261
41	Timer Circuit Board	1	SCEL0121	SCEL0121
42	Front label	1	CADS1563	CADS1563
43	Convex Head Philips M6 X 30mm	4	CPAR1313	CPAR1313
44	Rubber Foot Ø12.7 X 3.5mm	5	CBOR1249	CBOR1249
45	Noise Filter	1	SCEL0124	SCEL0124
46	Wire Connector (Wago)	1	CCNT1377	CCNT1377
47	Headless Philips Screw M2.5 X 20mm	4	CPAR1588	CPAR1588
48	O'ring	2	CBOR1262	CBOR1262
49	Headless Allen Screw M8 X 16mm	1	FPAR0713	FPAR0713
50	Spacer	4	FBCH1140	FBCH1140
51	Washer	4	FARR1141	FARR1141
52	Nut M2.5 DIN970	4	CPOR1629	CPOR1629

TECHNICAL SPECIFICATIONS

- Source: 127V or 220V AC 50/60 Hz
- Consumption: 750 W Average
- Heating Power: 450 W
- Motor Power: 1400 W
- Dimensions (WxDxH): 180x230x290 mm (without packaging), 200x250x300 mm (with packaging)
- Weight: 4.500 Kg (without packaging), 5.600 Kg (with packaging)
- Pollution Degree: 2
- Protection Degree: IPX1 (protected against spatter)
- Noise Level: Motor 70/90 dB @ 0.50M

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
Essix® SelectVac® Approximate Heating Times

Plastic	Thickness	Heating Time
Essix ACE® Plastic	.030"/¾mm	0:01:25
Essix ACE® Plastic	.035"/⅞mm	0:01:30
Essix ACE® Plastic	.040"/1mm	0:01:35
Essix® PLUS™ Plastic	.035"/⅞mm	0:02:20
Essix® PLUS™ Plastic	.040"/1mm	0:02:25
Essix C+® Plastic	.040"/1mm	0:02:30
Essix® Nightguard Laminate	.120"/3mm	0:02:30
Essix® Dual Laminate	.080"/2mm	0:02:45/0:00:50 (soft/hard)
Essix® Dual Laminate	.120"/3mm	0:03:00/0:01:00 (soft/hard)
Essix Tray Rite® Plastic	.020"/½mm	0:01:25
Essix® Model Duplication Material	.040"/1mm	0:01:30
Essix® Model Duplication Material	.060"/1½mm	0:01:35
Essix® Model Duplication Material	.080"/2mm	0:01:40
Essix A+® Plastic	.020"/½mm	0:01:20
Essix A+® Plastic	.030"/¾mm	0:01:30
Essix A+® Plastic	.040"/1mm	0:01:35
Essix A+® Plastic	.060"/1½mm	0:02:00
Essix A+® Plastic	.080"/2mm	0:02:05
Essix A+® Plastic	.120"/3mm	0:03:20
Essix® Sports Mouthguard Material	.120"/3mm	0:02:15
Essix® Sports Mouthguard Material	.160"/4mm	0:02:50
Essix® Laminated Sports Mouthguard Material	.160"/4mm	0:02:30
Essix® Laminated Sports Mouthguard Material	.200"/5mm	0:03:30

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Dentsply Sirona Orthodontics Inc.
Obelis S.A