# Hilliard Thermoplier® Pliers

DIRECTIONS FOR USE  
INSTRUCCIONES DE USO  
MODE D'EMPLOI  
ISTRUZIONI PER L'USO  
GEBRAUCHSANWEISUNG

Manufacturer:  
**Dentsply International Raintree Essix**  
7290 26th Court East  
Sarasota, FL 34243 USA  
Tel: 800.883.8733  
[www.essix.com](http://www.essix.com)

#MKHTP RTE-179-18 Rev.7 9/18

**R<sub>x</sub>**  
ONLY

# Hilliard Thermoplier® Pliers

ENG

## Function:

These pliers are used to create or modify force points, to create space, to create slots or to form hooks in plastic aligners.

**Contraindications:** None known.

## ⚠ Warnings:

Treatment range using the Hilliard Thermoplier® Pliers is not to exceed 1.0mm of movement per month per tooth. Excessive rate of movement may cause root resorption.

Use caution with the dental burner and heated pliers, as severe burns may occur if handled improperly.

## Precautions:

For best performance, the pliers must be heated to the proper temperature for the plastic being used. Too much heat will result in a thin, weak section of the plastic. Too little heat will not fully thermoform the plastic, causing it to move back to the original position over time.

As treatment progresses and cumulative force point sizes increase, a maximum cumulative force point should not exceed 3.0mm as the surrounding material may weaken and become less effective. All Hilliard Thermoplier® Pliers indicated for force points can be adjusted to form a force point of up to 3.0mm in depth with the exception of the Thermoplier® Micro-Ramp Pliers which can be adjusted up to 2.0mm.

**Adverse Reactions:** None known.

## Step by Step Instructions:

1. Pliers must be sterilized before each use. See sterilization guidelines below.
2. Using the provided hex screw, adjust the pliers to the desired depth.
3. Using the dental burner, heat the designated side of the pliers to the appropriate temperature depending on the Essix® Brand Plastic being used:

**Essix ACE® Plastic 245°F / 118°C      Essix C+® Plastic 200°F / 93°C**

4. Measure the temperature of the heated pliers using a thermometer.
5. Place the heated pliers in the desired position on the clear aligner and gently squeeze together.

## Thermoplier® Pliers Purpose/Indication:

Plier	Create Force Points	Purpose/Indication	Position on Aligner
Thermoplier® Undercut Enhancer Pliers (Red)	No	Adjust loose aligners by tightening/ extending plastic at interproximal points	Undercut areas/ embrasures on both sides of the aligner alternating undercuts
Thermoplier® Maxillary Tooth Moving Pliers (White)	Yes	Tip or rotate larger teeth	Center of the target tooth
Thermoplier® Mandibular Tooth Moving Pliers (Blue)	Yes	Tip or rotate smaller teeth	Center of the target tooth
Thermoplier® Micro-Ramp Pliers (Yellow)	Yes	Tip, torque or rotate smaller teeth laterally; Allows contact close to the incisal edge	On the correct location on the target tooth
Thermoplier® Mesial-Distal Pliers (Purple)	Yes	Lateral movement of both posterior and anterior teeth; Allows contact close to the incisal edge	Appropriate interproximal areas; a force point is placed on both sides of the aligner, angled to push the desired tooth mesial or distal
Thermoplier® Large Circle Pliers (Gray)	Yes	Increase the size (diameter) of a large existing force point; Elevate H.I.T.® Hooks off of teeth	In the area around the force point or hook
Thermoplier® Small Circle Pliers (Pink)	Yes	Increase the size (diameter) of small existing force point	Location of the existing force point
Thermoplier® Bite Plane Pliers (Black)	No	Make bite plane aligners without acrylic	On the resulting mark from articulating paper
Thermoplier® Bubble-Forming Pliers (Brown)	No	Create additional space in the aligner for the target tooth to move into	The location opposite the force point on the target tooth
Thermoplier® Hook-Forming Pliers (Green)	No	Produces the projection needed to make a hook for the attachment of elastics	Center of the target tooth
Thermoplier® Elastic Slot-Forming Pliers (Orange)	No	Forms the slot in the projection to hold elastics	On projection formed by the Thermoplier® Hook-Forming Pliers

Reference the MTM® Clear•Aligner In-Office Quick Start Guide for detailed examples of use.

## Cleaning and Sterilization Guidelines:

Clean, rinse and dry before sterilization.

Sterilize by autoclave/moist heat at 121°C / 250°F for 30 minutes (gravity displacement cycle). Minimum recommended drying time is 30 minutes.

## Maintenance Suggestions:

**AFTER** sterilizing, lubricate **SPARINGLY ONCE A WEEK** with dental handpiece oil. One drop or light spray is sufficient; wipe off excess lubricant. The major cause of staining or discoloring is **OVERUSE** of lubricants or “sticky joint” treatments. Lubricants, waxes and “sticky joint” remedies will vaporize inside the chamber and saturate all other items placed inside.

If you are experiencing brown or yellow staining, scrub with **ISOPROPYL** alcohol and steel wool (type A extra fine).

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