Chewies™ Aligner Tray Seater

DIRECTIONS FOR USE
INSTRUCCIONES DE USO
INSTRUÇÕES DE UTILIZAÇÃO
MODE D'EMPLOI
ISTRUZIONI PER L'USO
GEBRAUCHSANWEISUNG
GEBRUIKSAANWIJZING
BRUKSANVISNING
BRUGSANVISNING

Manufacturer:

Dentsply Sirona Orthodontics Inc.

7290 26th Court East

Sarasota, FL 34243 USA

Tel: 800.883.8733

EC

REP

Emergo Europe
Prinsessegracht 20
2514 AP, The Hague
The Netherlands

#MKCATS RTE-132-18 Rev.4 7/18

Chewies™ Aligner Tray Seater

ENG

Intended Use: Intended to be used to eliminate the air space between the thermoformed aligner/retainer and the teeth in order to properly seat the aligner/retainer.

Indications for Use: Indicated for patients wearing thermoformed aligners or retainers.

Contraindications: None known.

 **Warnings:** Use of excessive pressure to seat the aligner or retainer may cause damage to it.

Precautions: Momentary use for seating aligner or retainer only.

Seat only one aligner or retainer at a time. Trying to seat two at once may result in damage to the aligners or retainers.

Adverse Reactions: None known.

Step by Step Instructions:

1. Place the aligner or retainer over the teeth by hand and apply slight pressure until it is seated most of the way.
2. Insert the Chewies™ Aligner Tray Seater into your mouth and position it between the upper and lower arches.
3. Gently bite down on the Chewies™ Aligner Tray Seater and work the aligner or retainer into place until fully seated.

Safe Disposal of device:

Follow your country-specific laws, directives, standards and guidelines for the disposal of used devices.

Ensure that the parts are not contaminated upon disposal.

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